



# Colloidal Gold Nanospheres 20, 40, 60, 80nm

## APPLICABLE TO

GNS-20-ODI, GNS-20-OD5, GNS-20-OD10  
GNS-40-ODI, GNS-40-OD5, GNS-40-OD10  
GNS-60-ODI, GNS-60-OD5, GNS-60-OD10  
GNS-80-ODI, GNS-80-OD5, GNS-80-OD10

## SPECIFICATIONS

Particle Size	Diameter	% CV	% Odd Shapes
20nm	20nm±1nm	<8	<5
40nm	40nm±3nm	<8	<5
60nm	60nm±3nm	<8	<5
80nm	80nm±5nm		<10

## PRODUCTS

GNS-20-ODI, GNS-20-OD5, GNS-20-OD10  
GNS-40-ODI, GNS-40-OD5, GNS-40-OD10  
GNS-60-ODI, GNS-60-OD5, GNS-60-OD10  
GNS-80-ODI, GNS-80-OD5, GNS-80-OD10

## SAFETY

See MSDS

### Limited Use License

Use and Use restrictions. The products are sold, and deliverables of any services are provided for the purposes of the buyer's internal in vitro research, development or educational use only, not for in vivo, or any therapeutic or diagnose use, nor for resale, or for providing services or any other commercial use of any kind, including without limitation, for any trGNSfer in any form (including as part of a kit) to a third party, for analysis or reverse engineering of the Product or for manufacturing. Products should only be use in accordance with any safety data sheets, guidance or protocols that we issue from time to time and are available for download from our Site. Protective clothing should be used at all times when handling our products. Safety data -sheets relating to all Products are available for download from the Site or upon request.

Lateral Dx grants no other license or rights under any intellectual property in respect of Products to services deliverables and in particular grants no license to use nay Product or deliverables for any commercial purposes.

Sale of Products or service deliverables by us or our authorised distributors are expressly conditional upon the customer's agreements with these restrictions, which the customer gives upon placing an order for products or deliverables. If you wish to use any Products or deliverables for any purpose other than your own internal research as described above, you will require an additional license from Lateral Dx. Please contact [technical@lateraldx.com](mailto:technical@lateraldx.com)

## INTRODUCTION

Lateral Dx's colloidal gold nanospheres are high quality, with a narrow size distribution and uniform shape, supplied in a range of sizes and concentrations for maximum flexibility.

## STORAGE

The colloidal gold nanospheres are shipped at ambient temperature, upon receipt store between 2-8°C, do not freeze.

## DILUTION

Dilute in ultrapure water